



Perfect mixing quality with 30 % less wear!

We tested: Replica fodder mixing wagon knives compared to Strautmann original knives with hardened blade
Test setup: In a Verti-Mix 2401 Double, a mixing auger was equipped with replica knives.
 The second mixing auger was equipped with Strautmann ORIGINAL knives.
 After 2-year use with about 15,000 t of mixed fodder, we evaluated the result:


REPLICA

No more cutting performance!

Considerable abrasion

Serrated edges
hardly existing

Knife tip <1 mm
of residual thickness



- Frequent change of knives
- Poor mixing result
- Inadequate supply of herd


STRAUTMANN GENUINE PARTS

Cutting performance still existing!

Considerably less wear

Serrated edges
undamaged

Knife tip 3 mm
of residual thickness



- Cost-effective
- Consistent mixing quality
- Optimum supply of herd